

## **Maxpar Human Immune Activation IMC Panel Kit**

Catalog number: 201504

## Storage:

Antibodies: 4 °C. Do not freeze.
Cell-ID™ Intercalator-Ir: -20 °C

Reactivity: Human

**Application:** IMC paraffin **Reactivity:** Human

Target	Clone	Tag
Granzyme B	EPR20129-217	167Er
Ki-67	B56	168Er
PD-1	EPR4877(2)	165Ho
PD-L1	SP142	150Nd
Nucleic acid		191Ir/193Ir

## **Contents:**

- Cell-ID<sup>™</sup> Intercalator-Ir (125 μM; 500 μL)
- Maxpar antibodies\*, 25 μg, 0.5 mg/mL (See table for panel.)
- \* The antibodies are provided in individual tubes, not a pre-mixed cocktail.

## **Technical Information**

**Description:** The Maxpar® Human Immune Activation Panel IMC Kit enables 5-plex phenotyping of activated T cells in formalin-fixed, paraffin-embedded (FFPE) tissue sections using Imaging Mass Cytometry $^{\text{TM}}$  (IMC $^{\text{TM}}$ ). This panel kit includes 4 antibodies and a nucleic acid stain, and it leaves 31 open detection channels. This kit can be combined with the Maxpar Human Tissue Architecture and Maxpar Human Tumor-Infiltrating Lymphocytes IMC Panel Kits.

**Recommended usage:** For staining with the Maxpar Human Immune Activation IMC Panel Kit, prepare human FFPE tissue using standard techniques according to the Imaging Mass Cytometry Staining Protocol for FFPE Sections (PN 400322). Additional materials and equipment may be required for tissue staining and acquisition.

For technical support visit techsupport.fluidigm.com. | For general support visit fluidigm.com/support. For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legal/notices. **Limited Use Label License:** The purchase of this Fluidigm Instrument and/or Consumable product conveys to the purchaser the limited, nontransferable right to use with only Fluidigm Consumables and/or Instruments respectively except as approved in writing by Fluidigm. **Trademarks:** Fluidigm, the Fluidigm logo, Cell-ID, Imaging Mass Cytometry, IMC, and Maxpar are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2019 Fluidigm Corporation. All rights reserved. 08/2019